



EB Portable Chiropractic Table User Manual

OUR
BRANDS:



Applies to the following Models:
PCT1000



Important! This User Manual contains important information for the user of the product. All who use this product should review and completely understand its contents. Remember to keep this manual in a place where it is always available to those using the product.

Table of Contents:

Introduction:	2
Safety Precautions & Definitions:	2
Installation & Unpacking:.....	3
Technical Specifications:	4
Operation:	5
Preventative Maintenance:.....	6
Upholstery Cleaning & Care:	6
Service & Replacement Parts:.....	7
Replacement Part Information:	8
Warranty:	9
Manufacturer Information.....	10






Introduction:

The EB Portable Chiropractic Table is the perfect design for portability and function in a chiropractic table. It features a lightweight design with fixed arms and fixed height. The EB Portable Chiropractic Table is strong, stable, and easy to transport and set up. This table provides a viable alternative to bulky stationary designs. In its folded state, the PCT1000 is easily transported or stored in the smallest places.

Intended Use: The EB Portable Chiropractic Tables are intended for use in a clinical environment to support patients during treatment.

Safety Precautions & Definitions:

Warning and Precautions through this manual are indicated by specific symbols. Applicable symbols are shown below along with a description of the symbol. Review these symbols and all safety precautions before operating the table.

	CAUTION: Text with a “CAUTION” indicator will explain possible safety infractions that could have the potential to cause injury.
	WARNING: Text with “WARNING” indicator will explain possible safety infractions that will potentially cause serious injury and equipment damage.
	Pinch Point: Indicates space between manually adjustable moving and stationary parts where body parts may become caught, leading to minor injury.
	Instructions for Use: Indicates the need for the user to consult the instructions for use.
	Indicates the device manufacturer.

Before using this equipment follow the safety precautions below:



WARNING: Never leave patient on table unattended.



WARNING: Under no circumstances should the table be modified from its original design. Contact Pivotal Health Solutions if you need help with service.

- Before using this equipment make sure you have read and understand this entire manual.
- Explain proper use to all employees. This table should be used only by qualified, trained personnel.
- Never exceed the table's maximum supported weight.
- Use only on dry surfaces. Ensure the tables feet do not slide.
- This table is NOT for patient transport.
- NOT intended for use with accessories other than those indicated in this manual.
- Always explain to the patient how to properly get on and off the table. Be sure to assist patient to prevent falls.
- Keep patient's hands positioned on the top side of cushions at all times.
- Keep hands and feet away from moving parts and pinch points.
- Ensure controls have locked properly prior to positioning or releasing patient.
- Never leave patient unattended on table.
- Follow preventative maintenance instructions provided in this manual.
- Ensure all components have been inspected for damage and are fully functional.
- Inspect cables and connectors before each use.
- If you need service assistance contact Pivotal Health Solutions at 1-800-743-7738.

Installation & Unpacking:



WARNING: Do not install/use this equipment adjacent to, or stacked with, other equipment as it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify they are operating normally.

UNPACK AND INSPECT TABLE IMMEDIATELY UPON DELIVERY – DO NOT SIGN THE SHIPPING RECEIPT UNTIL YOU HAVE CONFIRMED THE PRODUCT IS FREE OF DAMAGE. FAILURE TO NOTE SHIPPING DAMAGE ON THE SHIPPING RECEIPT WILL VOID THE WARRANTY.

1. After removing the shipping box, you may notice metal bands and plastic cable ties. These are used to secure various sections of the table and must be removed.
2. Once the table is out of the shipping box place table in a level location.
3. Inspect the table cushions, frame, hardware, etc. for damage.
4. If no damage is noted, check that any pre-installed hardware is tighten. Tighten if needed.
5. Contact Pivotal Health Solutions or your dealer with any questions or concerns.

Technical Specifications:



WARNING: Keep table out of high moisture areas. For indoor use only.

Operating Conditions: This table should be operated, between 40° F (4.4° C) and 100° F (37.8° C), with relative humidity ranging from 30% - 80%.

Transport & Storage Conditions: This table should be transported and stored in temperatures between -40° C (-40° F) and 70° C (158° F) with relative humidity 10%-100%.

EB Portable Chiropractic Table (PCT1000)	
Overall Dimensions:	68" L x 22" W x 19" H
Folded Dimensions:	34" L x 9" W x 21.75" H
Table Weight:	35 lbs.
Weight Capacity:	350 lbs.
Height:	19" H (Fixed)
Foam:	2" High Density Foam

Expected Life: This product has an approximate expected life time of 10 years when correctly handled, serviced, and inspected in accordance with these instructions.

Design and Quality: Pivotal Health Solutions is certified in accordance with ISO 13485, the standard for the medical device industry.

Product Changes: Pivotal Health Solutions, Inc. reserves the right to make products changes without prior notice. Contact your Pivotal Health Solutions representative for advice and information about product upgrades.

Operation:

Table Set Up:

1. After removing table from packaging, place the table on a flat, clean surface, and unbuckle the buckle at the head end of the table.
2. Unfold the table with the bottom side facing up.
3. Extend all of the legs into an upright position.
4. Once legs are extended attach carabiner clips onto all three legs, clip into the hole that intersects both leg sections (See photo below).
5. Turn table in upright position and flip arm rests down to use.




Table Tear Down:

1. Place the table on a flat, clean surface with the cushion side facing down. Undo the carabiner clips and let the lanyard hang straight down.
2. Fold the middle leg down flat first, then fold the end on top of the middle leg followed by folding the leg at the head of the table.
3. Verify that the carabiner clip lanyard is placed underneath the leg and not pinched in the mounting bracket. (See photo below)
4. Fold the table in half and secure with the buckle. Place table on side in transport/storage position.

Carabiner clip lanyard placed underneath leg



Carabiner clip lanyard placed underneath leg

Assisting Patient onto Table:

1. The practitioner or trained staff should always assist the patient on the table. Do not allow the patient to get on the table by him/herself.
2. The patient should check clothing for keys, wallet, belt, or any other items that may cause discomfort during treatment or damage upholstery.
3. Never leave patient unattended while on the table.

Assisting Patient off of Table:

- 1) Caution the patient not to get up from the table alone so as to prevent any strain that could void the treatment, and to ensure the patient is not lightheaded from lying down which can cause imbalance.
- 2) Allow patient to get full balance before coming to a standing position. Best practice is to keep a hand on the patient at all times.

Preventative Maintenance:

INSPECTION:

At least monthly the table should be thoroughly inspected for wear and tear, loose hardware and parts, and other damage.

- Inspect the table to make sure that there is no external damage or loose hardware. Tighten loose hardware. Hardware that does not tighten needs replacement.
- Inspect table for obvious signs of damage or wear such as cracked welds, loose bolts, frayed or damaged cords. Do not use a table with damaged cords.
- Check moveable components to ensure functionality.
- Check components to ensure they lock and function properly.
- **Contact Pivotal Health Solutions with questions or concerns.**

Upholstery Cleaning & Care:



CAUTION: Read carefully, improper cleaning will void the warranty. Never use alcohol-based cleaning agents.

DISINFECTING: There is a disinfecting wipe on the market called *Protex Ultra Disinfectant Wipes*. Many of our customers use these.

FOR DAY TO DAY CLEANING: A solution of 10% mild household liquid dish soap with warm water, applied with soft damp cloth. Rinse with clean water and dry.

FOR STAINS: Dampen a soft white cloth in a one to one (1:1) solution of Fantastik® and water OR Formula 409® and water. Rub gently and rinse with a water dampened cloth.

FOR MORE DIFFICULT STAINS: Dampen a soft white cloth with a solution of household bleach (10% bleach / 90% water). Rub gently and rinse with a water dampened cloth to remove bleach concentration.

WHAT NOT TO USE: Using the wrong cleaning agents you will VOID YOUR WARRANTY and crack, dry out and destroy your vinyl. Do not use any cleaning agents that contain alcohol, harsh chemicals or abrasives.

Service & Replacement Parts:

Refer to Warranty section in this manual for specific warranty information.

When replacing worn parts, use genuine Pivotal Health Solutions, Inc. parts by contacting our service department. When ordering replacement parts:

- Ask for the service department, 800-743-7738
- Have the model number and serial number available
- Specify parts by the numbers/description below

Instructions for Damaged Merchandise: The contents of this shipment have been checked and packed by experienced personnel. If your company arranged the shipping, damage should be noted on the Bill of Lading before signing to accept the shipment, then articles damaged in transit should be reported at once to delivering carrier and a claim must be filed by you. If freight was arranged by Pivotal Health Solutions, note any damage on the Bill of Lading before signing acceptance of product and notify Pivotal Health Solutions right away. Failure to inspect for shipping damage prior to signing for the product may affect your warranty. We can neither be responsible for, nor can we accept the return of merchandise damaged in transit during shipment arranged by outside parties.

Technical Information:


Pivotal Health Solutions will make available, when needed, circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist service personnel.

Resale or Disposal:

A chiropractic table is considered a medical device by the Food and Drug Administration. Therefore, it is necessary that Pivotal Health Solutions, Inc. be notified if the table is sold, destroyed, or otherwise disposed of. Please notify in writing providing your name and the serial number of your table to:

Pivotal Health Solutions
Attn: Quality
3003 9th Ave. SW
Watertown, SD 57201

Replacement Part Information:



EB Portable Chiropractic Table	
<u>Part Name</u>	<u>Part Number</u>
1" Nylon Web Handle	A2216
Coiled Lanyard w/Carabiner	TT0132
T Bench Legs	TT0012
1" Rubber End Cap	TT0048
1" Squeeze Release Plastic Buckle	TT0011
Face/Chest Cushion	TT0003-xxx
Pelvic/Leg Cushion	TT0006-xxx
Arm Cushion	TT0009-xxx

*xxx is a place holder for the color code. Contact Pivotal Health for the appropriate color code for your table.

Warranty:

Pivotal Health Solutions warrants to the **original purchaser** that each PHS Chiropractic product shall be free from defects in material and workmanship and agrees to remedy any such defect or to furnish a new or equal part in exchange therefore. Warranty covers components and repairs (labor) indicated below, commencing on date of original shipment. Labor coverage may not include replacement or installation of small parts or components. Pivotal Health Solutions Service and/or Sales will determine the technical level of service needed and approve labor coverage on a case by case basis.

Warranty Period: 1-year structural frame, electrical components, foam and upholstery, moving components, hardware, and accessories. Parts Only Warranty.

Pivotal Health Solutions reserves the right to inspect claimed defective products. No returns, replacements or repairs may be made without prior written consent. Please do not return any product to Pivotal Health Solutions or its service center without first contacting Customer Service for a Return Authorization number. Products returned without a Return Authorization number will be refused back to the customer at their expense.

Return Authorization can be obtained by contacting the service department at Pivotal Health Solutions. Refer to Service section of this manual for contact information.

This Warranty covers defects in materials and workmanship provided the product is used for the normal purpose for which intended. The customer shall be responsible for the proper set-up and use of the product. This Warranty does not include damage resulting from accident and expressly excludes normal wearing of parts or deterioration due to normal wear and tear, damage caused by improper set-up or storage, defect caused by transportation, accident, fire, flood, alteration, abuse, misuse, or negligence. Pivotal Health Solutions shall not be liable for any direct or indirect, incidental, exemplary, or consequential damages or delay, including damages for loss of income or loss of use.

Any obligation of Pivotal Health Solutions under this Warranty shall automatically and immediately terminate, without notice from or any further action by Pivotal Health Solutions and Pivotal Health Solutions shall have no responsibility for damages of any kind as a result of the occurrence of any of the following:

- accident, misuse, abuse or negligent use of the Products or any component thereof;
- any repair or alteration of the Products or any component thereof made outside Pivotal Health Solutions' authorization;
- improper installation, setup, or operation (including both mechanical and electrical) of the Products or any component thereof, which includes the operation of the Product not in accordance with the Product's operating manual;
- failure to provide normal maintenance for the Products or any component thereof in accordance with the Product Operating Manual.
- Alteration or obliteration of any identifying marks.

Note: Warranty periods listed apply only to products purchased in new condition. Any tables sold as a factory blemish or in used condition are sold 'as-is' and are limited to a (1) year warranty on the frame only, commencing on the date of delivery.

Manufacturer Information



Pivotal Health Solutions, Inc.
3003 9th Ave. SW
Watertown, SD 57201
USA



(in US) 800-743-7738



(international) 011-1-605-753-0110



Fax: 605-882-8398



Email: service@pivotalhealthsolutions.com



World Wide Web: www.pivotalhealthsolutions.com



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Rev A 04/20/2020